GMOD Conference Call Minutes Addendum 3pm April 14, 2016

- 1. Test stand re-qualification
 - Labs are required to complete 1 test on reference oil assigned by TMC to re-qualify stand used in precision matrix.
 - Y_i alone will be used for re-qualification. This test will be the first chartable test. R_i, Z_i, and Q_i will begin control charting with the first reference test after requalification.
- 2. Test Monitoring System updated
 - Lower control limits for R_i, and Q_i have been deleted. Only upper control limits will impact a stand's status for these 2 parameters.
 - Document can be accessed at https://dexostrc.org/reference. Select Information Letters and GMOD under Choose Filter. Select GMOD information Files. Since the test monitoring system has not yet started, the latest version is still called "Issue 1".
- 3. New WPD values
 - WPD values from the precision matrix calculated using the new weighting factors are available at http://ftp.dexostrc.org/reference/docs/GMOD/MatrixStatusReports/
- 4. Post-test ring end gap measurement
 - Post-test top and second ring end gap measurements will not be required.

Submitted by Andrew Buczynsky, 4/22/16